



VEAL INFUSION

Cat# 6009

<u>TEST</u>	<u>SPECIFICATION</u>
Appearance	Dry, fine, homogeneous, light to medium tan brown, Free of extraneous material
2% Solution	Light to medium dark tan yellow
Color (% Trans @450 nm)	15% or greater
Solubility (Held 5 days)	Clear to slightly hazy with no more Than slight trace sediment
pH	6.5 to 7.5 at 25 ⁰ C @
Total Nitrogen (Kjeldahl)	10% or greater
Ash (Residue on Ignition)	15% or less
Loss on drying	6% or less
Microbiological Analysis	
Growth on	
Peptone agar	Satisfactory
HI Broth	Satisfactory
Total plate count:	<5000 col/g
Yeast & Molds:	<100 col/g
Coliforms:	Negative
Salmonella:	Negative